## **Technical note**



## Apogee Driver Installation and verification test

### Applicable to Apogee Camera Range (Alta/Aspen/Ascent)

This guide will demonstrate how to install and test the Apogee camera drivers.

#### Initial considerations

For demonstration and verification purposes, this document has been created using an Aspen camera on the *Apogee Camera Control Demo Program*, freely available from the Andor website <sup>[1]</sup>. The following process will be similar for other applications.

#### Method

- Download the Apogee SDK from the Apogee downloads section of the Andor website <sup>[2]</sup>. The Apogee SDK is 32-bit but includes both 32-bit and 64-bit versions of the Apogee camera drivers.
- Extract the Apogee SDK .zip folder and run the setupApogeeSoftware.exe.
- Follow the on-screen instructions to install the software.
- Power ON the camera and connect to the PC using a USB cable. The camera should automatically be detected and the driver installed.
- Run the Apogee Camera Control Demo application to confirm camera operation <sup>[1]</sup>.
- Select Discover from the Camera menu to open the Camera and Filter Wheel Selection Window.

(Note, if setting up an Ethemet connection, see the Simple setup of Ethernet control for Apogee Aspen cameras document on the Andor Support webpage <sup>[3]</sup>)

🖳 Аро	gee Camera	a Control De	emo			
File	Options	Camera	Filter Wheel	View	Help	
1		Con	nect			
		Disco	over		Y	 _
lma	iging Mc	Disco	onnect		Γ	
	Mode Si	IO Po	ort Assignments	5		
Sequence	ce Count		1		·	

## **Technical note**



- Select USB 2.0 as the Connection Interface and then press the Search button. Any device that is detected will be listed in the Search for Device window. Select the device you wish to control and press OK.

Connection Interface	F Ethernet	Broadcast 172	26 . 39 . 255
Search for Devices			
Name	Interface	Id One	Id Two
	ush		
Aspen-230			
Search			

- If the connection is successful then the camera model will be listed in the *Camera Info* section in the bottom left corner of the program Window. The CCD temperature will also be displayed indicating that the drivers have successfully been installed and the program can successfully communicate with the camera

AD Gain 1 47 AD Offset 1 300 AD Gain 2 47 AD Offset 2 285 Camera Info S Driver Version: 55.32.3218 Sensor Size: 2048 x 2048 Fernivare Version: 306 Camera Model: Aspen-230 USB FW Version: a11 Filter Wheel Control S Model Model Postion	P Apogee Camera Control D	smo	
Binning X 1 + Binning Y 1 + Binning Y 1 + Readout Control AD Gan 1 47 + AD Gfiset 1 300 + AD Gan 2 47 + AD Gfiset 2 226 + Camera Info Driver Version: 55.32.3218 Sensor Size: 2048 x 2048 Ferrivare Version: 306 Camera Model: Aepen-230 USB FW Version: all	File Options Camera	Filter Wheel View Help	
Binning X 1 + Binning Y 1 + Binning Y 1 + Readout Control AD Gan 1 47 + AD Gfiset 1 300 + AD Gan 2 47 + AD Gfiset 2 226 + Camera Info Driver Version: 55.32.3218 Sensor Size: 2048 x 2048 Ferrivare Version: 306 Camera Model: Aepen-230 USB FW Version: all	🗃 🗄 🖳 🙆	🗣 (i) 🖬 🔿	
Binning Y 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	C-0 (0.00)		
Readout Control   peline Download   Off   AD Gain 1   47 (±)   AD Offset 1   300 ÷   AD Offset 2   285 ÷   Camera Info   Sensor Size:   2048 x 2048   Femware Vension:   55.32 3218   Sensor Size:   2048 x 2048   Femware Vension:   306   Camera Model   Appen-230   USB FW Vension:   a11		1000	
peline Download Off Oni AD Gain 1 47 AD Offset 1 300 AD Gain 2 47 AD Offset 2 285 Camera Info S Driver Version: 55.32.3218 Sensor Size: 2048 x 2048 Femware Version: 306 Camera Model Agen-230 USB FW Version: a11	Deriving 1		
AD Gain 1 47 AD Offset 1 300 AD Gain 2 47 AD Offset 2 285 Camera Info S Driver Version: 55.32.3218 Sensor Size: 2048 x 2048 Fernivare Version: 306 Camera Model: Aspen-230 USB FW Version: a11 Filter Wheel Control S Model Model Postion	Readout Control		
AD Offset 1 300 A AD Gan 2 47 A AD Offset 2 285 A Driver Version: 5.5.32.3218 Sensor Size: 2048 x 2048 Fernivare Version: 306 Camera Model: Appen-230 USB FW Version: a11 Filter Wheel Control A	peline Download	06	
AD Gan 2 47 2 AD Offset 2 285 2 Camera Info Driver Vension: 55.32.3218 Sensor Size: 2048 x 2048 Ferniware Vension: 306 Camera Model: Appen-230 USB FW Version: a11 Filter Wheel Control	AD Gain 1	47 💠	
AD Offset 2 285 C Camera Info S Driver Version: 55.32.3218 Sensor Size: 2048 x 2048 Finnware Version: 306 Camera Model Appen-230 USB FW Version: a11	AD Offset 1	300 ≑	
Camera Info Driver Vension: 5.5.32.3218 Sensor Size: 2048 x 2048 Femware Vension: 306 Camera Model Appen-230 USB FW Vension: a11 Filter Wheel Control Model Postion	AD Gan 2	47 😂	
Driver Version: 55.32.3218 Sensor Size: 2048 x 2048 Firmware Version: 306 Camera Model: Aspen: 230 USB FW Version: a11 E Filter Wheel Control S Model Model Postion Version Versio Version Version Version Ver	AD Offset 2	286 0	
Sensor Size: 2048 x 2048 Firmware Vension: 306 Camera Model: Aspen:230 USB FW Version: a11 Filter Wheel Control	Camera Info		
Firmware Version: 306 Camera Model: Aspen-230 USB FW Version: a11 E Filter Wheel Control R Model Model Postion Postion	Driver Version: 5.5.32.3218		
Camera Model: Aspen-230 USB FW Version: a11 E Filter Wheel Control F Model: Model: Postion Postion F	Sensor Size: 2048 x 2048		
USB FW Version: a11 I Filter Wheel Control Postion Postion I I I I I I I I I I I I I I I I I I I	Firmware Version: 306		
Filter Wheel Control			
Model Model Postion	USB FW Version: a11		
Postori Contra C	Filter Wheel Contro		
	Model Model		
	Position		
CD temp: 23.95 *C	CCD temp: 23.95 *C		

# **Technical note**



Further testing can be performed to check the camera functionality by taking an acquisition - set an exposure time and select the Take Image or Live Video buttons.

Apogee Camera Contro File Options Comm	Demo	
and the second sec	) <b>()</b> () () ()	
Imaging Mode C	ontrol 🛛	
Mode Single Sho	t •	
Sequence Count	1	
Sequence Delay	1.000 🕀	
Flename		
Exposure Contro	a 🛛	
Expsoure Type Light	-	
Exposure Time	1.00000 -	
Digitization Speed Normal	*	
Preflash Of	On	
Cooler Control		

## Troubleshooting

#### Camera not found

Confirm that the Apogee USB driver has been correctly installed from the Device Manager (Control Panel – System – Hardware – Device Manager)

Error message during installation process: 'ApogeeUSB.64bit mandatory prerequisite was not correctly installed'



Cancel the installation and disconnect all USB device (USB memory sticks, USB hard drives, etc - excluding keyboard, mouse). Run the Apogee SDK installer once again and complete the installation process ensuring that the Apogee camera has been attached to the PC using the USB cable. Any other devices can be reconnected once the installation has complete.

If you require further information or assistance, please contact your local Andor Product Support representatives <sup>[3]</sup>.

#### Useful links

<sup>[1]</sup> Apogee Camera Demo Program - <u>http://www.andor.com/pdfs/software/ApogeeCameraControlDemo.exe</u>

<sup>[2]</sup> Apogee SDK Download link - <u>http://www.andor.com/pdfs/software/apgSwInstall-x86.zip</u>
 <sup>[3]</sup> Product Support - <u>http://www.andor.com/ContactSupport</u>