

Simple setup of Ethernet control for Apogee Aspen cameras

Applicable to Apogee Aspen

This guide demonstrates how to setup Ethernet control of an Apogee Aspen camera using a network switch.

Initial considerations

If not already done so, download and install the Apogee SDK from the Andor website^[1]. Full instructions for installation of Apogee drivers can be found in the **Apogee Driver Installation and Verification test** document in the Support section of the Andor website^[2].

For demonstration and verification purposes, this document has been created using an Aspen camera on the *Apogee Camera Control Demo Program*, freely available from the Andor website ^[3]. The following process will be similar for other applications.

Ethernet Communication Setup

- Ensure the USB cable is disconnected from the camera.
- Using a network switch and Ethernet cables connect the switch, camera and PC together as shown in the diagram below before powering up the camera.



- Run the *Apogee Camera Control Demo Program*. The operating system may generate a security warning when starting the program. Select *Run* to allow the program to start.
- When the program starts, select Discover from the Camera menu to open the Camera and Filter Wheel Selection Window

🖳 Apoge	e Camera	a Control De	emo			_
File O	ptions	Camera	Filter Wheel	View	Help	
Tack I	-	Conr	nect			
		Discover		<u> </u>		
Imagi	ng Mo	Disco	onnect			
Mode Si		IO Port Assignments				
Sequence (Count		1			



- Select *Ethernet* as the Connection Interface and then press the *Search* button.

Any device that is detected will be listed in the *Search for Device* window. Select the device you wish to control and then press OK. **Be sure to take a note of the camera IP address listed in the** *Id One* field and in the selected device field.

USB 2.0	Ethernet	Broadcast Address: 172	. 26 . 39 . 255
Search for Devices			
Nama		140/10	Id Time In
Aspen-230	ethemet	172.26.35.73	80
Search	Vetwork Search		
Search	Vetwork Search		

If the connection is successful then the camera model will be listed in the *Camera Info* section in the bottom left corner of the program Window. The CCD temperature will also be displayed indicating that the drivers have successfully been installed and the program can successfully communicate with the camera.

程 Apogee Camera Control Demo	(()
File Options Camera Filter Wheel View Help	
🗃 🗏 🕼 🙆 💭 🕦 🖬 🔿	
Brining X 1 🚓	
Binning Y 1 💠	
Readout Control 🔯	
Ppeline Download Off On	
AD Gain 1 47 💠	
AD Offset 1 300 ≑	
AD Gan 2 47 🚭	
AD Offset 2 286 🗇	
Camera Info	
Driver Version: 5.5.32.3218	
Sensor Size: 2048 x 2048	
Firmware Version: 306	
Camera Model Aspen-230	
USB FW Version: a11	
Filter Wheel Control	
Model Model	



The camera IP address can be used to enable Ethernet control of the camera when using 3rd party software packages such as Maxim-DL.

If you require further information or assistance, please contact your local Andor Product Support representatives ^[2].

Useful links

^[1] Andor driver pack for Apogee cameras - <u>http://www.andor.com/pdfs/software/apgSwInstall-x86.zip</u> ^[2] Product Support - <u>http://www.andor.com/ContactSupport</u>

^[3] Apogee Camera Control Demo Program - <u>http://www.andor.com/pdfs/software/ApogeeCameraControlDemo.exe</u>.